

SpotOptics s.r.l. – leaders in accurate metrology

OMI-CMOS-DUV to NIR

VERSATILE WAVEFRONT SENSOR

- **Accurate metrology in single pass**
- **Optical elements, lasers and laser diodes**
- **Test any focal length and diameter (with accessories)**
- **Large dynamic range**
- **For R&D and production**
- **Optimized for DUV->NIR wavelength ranges (170-1000nm)**



TECHNICAL SPECIFICATIONS

HARDWARE

Test	Optical elements, lasers and laser diodes
Power of laser diode that can be tested	Few mW. Higher powers require a power reduction system (available)
No of spots (see cameras below)	35x35 (DUV-VIS-NIR) for a pupil size of 7mm
Diameter and focal length of standard lenslet arrays	<ul style="list-style-type: none"> • $\phi=0.2\text{mm}$, $f=22\text{mm}$ – for UV-Vis region from 170-750nm • $\phi=0.2\text{mm}$, $f=11\text{mm}$ - for NIR region – from 751-1000nm • Other lenslet arrays with smaller diameter and/or shorter focal length are available on request

SOFTWARE

Software (control and analysis)	Sensoft for 64bit Win 10
--	--------------------------

ACCURACY AND REPETABILITY

Resolution	1nm
RMS repeatability of Zernike coefficients (zonal wavefront)	<2nm rms
RMS repeatability of modal wavefront measurements	< $\lambda/100$
Accuracy and dynamic range	$\lambda/20$ - $\lambda/100$ (calibration dependent), $\pm 50 \lambda$

CAMERA


Detector, wavelength range and cooling	CMOS BSI technology (DUV-VIS-NIR). Quantum efficiency: Over 40 % at 200 nm, < 70 % at 650 nm. Uncooled
Chip size, pixel size	1280 x 1024 pixels, each of 7.4 μm . 9.5 x 7.6 mm ²
Connection, A/D convertor bits	Gigabit Ethernet, 12-bits
Acquisition speed and Frame transfer	28 Hz, global shutter
Triggering	Yes: external Trigger, SW Trigger
Exposure time (min, max)	280 μs , 559 msec

ACCESSORIES

Light sources, beam expanders and compressors	High quality LD with lens at test wavelength, beam expanders/compressors
--	--

OMI DUV-VIS-NIR (from 193nm-1100nm)

Quantum Efficiency Image Sensor






Quantum efficiency curve of CMOS

Other details

- Resolution: 1280 x 1024 pixels
- Pixel size: 7.4 μm x 7.4 μm
- Chip size: 9.47 mm x 7.58mm
- Image rate: 28hz (full resolution)
- Max. exp. time: 559 msec
- Connection: 10GigaEthernet

Other details

- Resolution: 35x35 spots (max)
- Lenslet pitch and focal length (UV and VIS):0.2mm, 22mm
- Lenslet pitch and focal length (NIR):0.2mm, 11mm
- Calibration unit for parallel light: Static or motorized high-quality collimator with LD/LED at test wavelength
- Motor step: 2.5 μm



SENSOFT: THE SOFTWARE


Sensoft: The modular software package

- Fully controls the hardware of OMI
- Performs the Shack-Hartmann (SH) analysis
- Computes Zernike coefficients, diagnostics (alignment and correct focal plane), zonal and modal wavefront, MTF, spot diagram
- Has a Loop mode for on-line adjustment of optical systems

OMI in your production line:


- OMI – with its own PC - can easily be adapted to the production line
- It can work in a closed-loop with the PC of the manufacturing machine
- A software module defines the IP communication protocol and transfers the results between the PCs in the Local Area Network

ON-LINE ALIGNMENT IN A FAST LOOP



- The alignment of complex optical systems becomes easy by monitoring coma and astigmatism in a continuous loop
- The individual (x, y) components of coma and astigmatism, as well as the total coefficients are displayed
- The optimization can be done for one component at a time, as the software can display one component of interest

MTF MEASUREMENTS



MTF after subtracting the contributions of tilt and defocus present in the data.

OMI with CMOS-DUV CAMERA



OMI

- **f_l=22mm, ϕ = 0.2mm for UV-Vis (170-750nm)**
- **f_l=11mm, ϕ = 0.2mm for NIR (750-1000nm)**
- **Max. Resolution 35x35 spots for a pupil size of 7mm**

PHYSICAL

Camera:

CMOS, 10GigE (170-1000nm), Global shutter

Chip size: 1280x1024x7.4 μ pixels, 9.47.x7.58 mm²

Dimensions: 59 (L) x 59 (W) x 102 (H) mm

Weight: 460 gm

KEY FEATURES

Measurement technique

Shack-Hartmann wavefront sensor

Test in parallel light or at the lens focus in single pass

Parallel light (with a calibration unit)

At the focus of the lens (with pinhole calibration unit)

Light sources with different wavelength available

Calibration units available

High-quality parallel light source (motorized or manual)

Pinhole calibration unit

Accessories

Light sources and beam expanders/compressors. Collimators

SOFTWARE

- Full waterfront analysis: Zernikes, zonal and modal WF, Spot diagram, MTF, EE, PSF, M²
- Easy alignment of lens group via software: graphical indication for correction using coma and astigmatism
- Stabilization of lasers: graphical indication of focusing of laser beam